

# IN-GEL DIGESTION PROTOCOL

## *Solutions Used During Protocol*

25mM Ammonium Bicarbonate  
25mM Ammonium Bicarbonate / 70% ACN  
10mM DTT in 25mM Ammonium Bicarbonate  
50mM Iodoacetamide in 25mM Ammonium Bicarbonate  
10ng/ $\mu$ l Trypsin in 25mM Ammonium Bicarbonate  
50% ACN / 5% Formic Acid (or TFA)

Keratin contamination is a major hazard of in-gel digestion. Therefore, 'powder free' gloves should be worn at all times during the protocol, and caution should be taken at all stages to minimize risk of contamination; i.e. lids on tubes, tip boxes and solutions should be kept on as much as possible.

In order to minimize sample loss, siliconized (hydrophobic coated) tubes are used throughout the process.

## *Wash/Dehydrate Gel Pieces – to wash gel pieces and remove coomassie stain.*

1. Add 50 $\mu$ l 25mM Ammonium Bicarbonate /70% ACN to gel pieces and incubate for 10 min.
2. Remove supernatant using a gel-loading pipette tip.
3. Repeat steps 1 and 2

## *Reduction and Alkylation – to reduce disulphide bonds and modify cysteines to prevent bonds from re-forming.*

4. Add 30 $\mu$ l 10mM DTT in 25mM Ammonium Bicarbonate to dried gel pieces and incubate at 50°C for 45 minutes.
5. Remove supernatant using a gel-loading pipette tip. Add 30 $\mu$ l 50mM Iodoacetamide in 25mM Ammonium Bicarbonate and incubate in dark for 1hr.
6. Remove supernatant using a gel-loading pipette tip. Add 30 $\mu$ l 25mM Ammonium Bicarbonate / 70% ACN and incubate for 5 min, then repeat this step.

## *Digestion*

7. Add enough trypsin solution to barely cover gel pieces. For a 2D spot this will be about 5 $\mu$ l, for a 1D band 10-15 $\mu$ l. Allow gel pieces to swell for 5 minutes.
8. Overlay gel pieces with a further 30 $\mu$ l 25mM Ammonium Bicarbonate and incubate at 37°C for 4hrs (or overnight).

### *Peptide Extraction*

9. Remove and collect supernatant in a new siliconized tube using a gel-loading pipette tip
10. Add 30 $\mu$ l 50%ACN / 5% Formic Acid and incubate for 10 minutes.
11. Remove and collect supernatant using gel-loading pipette tip in same siliconized tube as step 11.
12. Speed-vac collected supernatant down to 5-10 $\mu$ l.
13. Add 5 $\mu$ l 0.1% Formic Acid.